



AQ400

DISCRETE ANALYZER FOR ENVIRONMENTAL TESTING

AQ400 is a flexible analyzer that uses the principle of discrete analysis where each test occurs in a separate or discrete reaction vessel.



AQ400

AQ400 is ideal when many and varied tests are needed on different samples and/or individual results are needed immediately.

USEPA, ASTM, ISO

and other international regulatory compliant methods are available.

METHODS INCLUDE

Alkalinity
Ammonia
Chloride
Cyanides
Nitrate/Nitrite
Nitrite
Phenol
Phosphate, ortho
Phosphorus, total
Silicate
Sulfate
Total Kjeldahl Nitrogen

PLUS MANY MORE



HOW DOES THE AQ400 WORK?

The AQ400 robotic sampling arm works in conjunction with a stepper motor-driven syringe that is responsible for aspirating, dispensing and mixing accurate and precise quantities of sample and reagent in miniaturized test tubes called reaction wells.

The sample and reagent mixture are incubated in heated reaction wells until the reaction is complete. A single aliquot is then transferred into a 10 mm path length optical glass cuvette where the absorbance is read. Each sample is read in the same cuvette, in the same position in front of the detector. This is similar to SEAL colorimetric flow systems, known for their high reproducibility and lowest detection limits.

The flow through cuvette eliminates the issue of reaction well variability and scratching found in direct read discrete systems. In the SEAL AQ series, liquid is moved not the cuvette – fewer moving parts equal higher stability and reliability of the system.

Once the absorbance is read, the glass cuvette is thoroughly cleaned and checked, ensuring no carryover or cross contamination.

▶ **As a market leader**, SEAL has over 1,000 applications available and under continual development. Markets include water, wastewater, soil, plant, fertilizer, food and beverage. Please contact us for your specific application.

AQ400 – The latest in discrete analyzer technology

Designed specifically for the environmental market

ADVANTAGES AND BENEFITS



- ~ 100% optical glass stop-flow cuvette
- ~ 10 mm optimum path length
- ~ Longer path lengths available




- ~ Low cost, disposable wells used for each discrete reaction
- ~ Constant heating and programmable reaction time ensure reaction reaches completion
- ~ Automatic reagent level sensing verifies sufficient reagent volume



- ~ Optional integrated cadmium coil reduction for nitrate/nitrite determination
- ~ Cadmium coil is sealed and valve controlled
- ~ In situ coil regeneration is fully automated



- ~ Unique probe washer for cleaning sample probe between sample and reagent
- ~ Eliminates cross contamination
- ~ Ideal for wastewater




- ~ Highly flexible software designed with user input
- ~ QCPro™ Data Quality System – allows the user to specify QC types, limits and corrective actions

- ▶ **True unattended operation – including ability to run overnight**
- ▶ **Automated standard preparation and dilution of over range samples**
- ▶ **Tests programmable per sample to reduce analysis time**
- ▶ **Add samples after a run has started**
- ▶ **Total volume per test only 500 – 600 µL**
- ▶ **Different size sample trays are available to accommodate different workloads**
- ▶ **Segregated chemical waste and wash minimizes waste disposal**
- ▶ **LIMS compatible – export in .csv format**



SIMPLE USER MAINTENANCE



The AQ400 can be used as a standalone spectrophotometer.

The optional vial adaptor can accept vials of different sizes for reading tests such as COD.

SEAL
Analytical



www.seal-analytical.com

Colorimetric Nutrient Analyzers

DISCRETE ANALYZERS



AQ270



AQ300



AQ400

SEGMENTED FLOW ANALYZERS



AA100



AA500



QuAAtro39

50 Years' Experience in Environmental Analysis Built into Every Analyzer

50 years' experience in designing, developing and manufacturing automated wet chemistry analyzers specifically for very low detection levels in environmental applications has helped SEAL to apply the most useful, easy to use features into the SEAL range of Discrete and Segmented Flow analyzers. The SEAL analyzers are widely acknowledged as the best for environmental analysis, giving you everything you need to achieve equal or superior results to the manual and approved laboratory methods the SEAL analyzer replaces.

Digestion Systems

FOR METALS AND TKN,
TP DIGESTION



BD50



SmartBlock II



DEENA 3

SEAL
Analytical



www.seal-analytical.com

SEAL Analytical is a global company with offices worldwide - contact us at:

SEAL Analytical, Inc.
6501 West Donges Bay Road
Mequon, WI 53092
United States
Tel: +1 (262) 241 7900
Fax: +1 (262) 241 7970
sales@seal-us.com

SEAL Analytical Ltd.
3 Talisman Business Centre
Duncan Road
Park Gate
Southampton, SO31 7GA
United Kingdom
Tel: +44 (0) 1489 864 400
sales.uk@seal-analytical.com

SEAL Analytical Netherlands
ROHASYS BV
Provincienbaan 4
5121 DL Rijen
The Netherlands
Tel: +31 161 240152
Fax: +31 161 240153
info@rohasys.com

SEAL Analytical GmbH
Werkstrasse 5
D-22844 Norderstedt
Germany
Tel: +49 (0)40 60 9292 9-00
Fax: +49 (0)40 60 9292 9-02
info.germany@seal-analytical.com

SEAL Analytical Shanghai
Room 413, 12th Building,
No. 128 Xiangyin Road,
Shanghai, 200433
China
Tel: +86 21 3362 5002
Fax: +86 21 3362 5002

SEAL Analyzers are monitoring environmental samples in every corner of the globe. They are manufactured in the USA, Germany and the Netherlands. Engineering and chemistry support is provided from SEAL global facilities in USA, Germany, England, the Netherlands and China along with a worldwide network of specialist distributors.

COMPREHENSIVE SUPPORT

We offer comprehensive applications, technical service and software support.

INCLUDING

- ▶ A choice of preventative maintenance and service contracts to meet your specific requirements
- ▶ In-house and online training
- ▶ Guaranteed availability of genuine consumables and spare parts
- ▶ Adaptation of methods to specific requirements such as matrix, range or detection limit
- ▶ Continuous in-house development of software to incorporate new customer requested features

Robotic Handling Systems

SEAL Robotic MiniLab systems for automating sample pretreatment in the laboratory — improving your sample handling efficiency. Typical applications include BOD, pH, COD, Alkalinity, and conductivity measurements with options such as decapping/capping, sample splitting, and filtration. Call us about your laboratory needs and we will design a robot to suit you.



SEAL MiniLab