

# VASCO



**Extremely wide sample concentration**

ppm orders

← **range capability** →

Several 10 %



*When knowing particle size counts !*

## IDEAL FOR

*Formulation stability  
Nanoparticle aggregation  
Emulsions dispersion  
Pharmaceuticals  
Petrochemicals  
Polymers  
Liposomes and bio-colloids  
Pigments and inks ... and more*

**CORDOUAN**

T e c h n o l o g i e s





## Dedicated to nanoparticle size analysis

The **VASCO™ particle size analyzer** is a unique instrument for nanoparticle suspension and colloidal characterization, based on enhanced Dynamic Light Scattering (DLS) technology. Thanks to a patented technology developed in collaboration with the French Institute of Petroleum (IFP), **VASCO™** is the most efficient solution for concentrated and/or opaque suspension characterization.



## Main characteristics

- Based on enhanced Dynamic light Scattering (DLS)
- Embedded sample cell with patented DTC system
- Particle size (diameter) : 0,5nm - 10 µm ;
- Sample concentrations : 0.1ppm - 40%w/wt ;
- NanoQ proprietary software dedicated to particle size
- On-line sample cell option → size kinetics study ;
- Optical filter option to improve measurements on fluorescent samples ;



## Technologies & innovations

VASCO™ : a **unique** particle size analyzer

- An embedded sample cell made of a silica prism ;
- Patented Dual Thickness Controller system for accurate sample thickness control ;
- Measurement of **dark / concentrated** samples without dilution.

## Innovative Cell design : simple, robust, artifact free


**Simple** : The DTC sample cell design allows simple and easy **sample preparation** and prevents any strong dilution. It is **compatible with organic solvent** and can measure **micro-volume** of sample. An upgradable option for **on-line measurement** can complete **study of kinetics**.

**Artifact free**: cell prevents A simple adjustment of **DTC** is sufficient to reduce the sample layer from 2 mm to 200 µm (ultra-thin layer sample's volume). This Dual Thickness Controller for reducing the measurement volume prevents problems due to multiple-scattering as well as local heating and guarantees reliable measurements of dark media or high particle concentrations. In comparison with cuvette cell concept, sample needs to be strongly diluted and by consequence influence the properties of representation on the size distribution !



# VASCO technology led to its best

Think « out of the box »



VASCO offers the extremely wide sample concentration range capability, from very dilute sample or poor contrasted particle to high concentrated and/or opaque formulation. Cordouan technologies unique optical design patented makes exceptionally wide sample concentration range measurements practical. This enables use in a wide range of applications an industries that needs to be close on the end product without dilution.

## Key benefits

- **Higher detection efficiency** in **opaque/dark media** , as well as in **diluted** ;
- Limited sample preparation
- Solvent-proof embedded sample cell : **no consumable** ;
- Proprietary Padé Laplace algorithm for polymodal sample analysis.
- Extended range concentration samples that limited time preparation.
- 2 to 3 decades concentration range improvement over classical DLS.



## High performances for advanced applications



Pharmaceutical



Cosmetics



Chemistry



Advanced colloids



Polymers

- NP synthesis and functionalization study
- Drug delivery optimization
- Quality control in manufacturing process
- Fundamental study of electrophoretic physics
- Cosmetic and industrial emulsion stability study
- Nanoparticle formulation and synthesis optimization
- Advanced colloidal stability analysis and optimization
- Ink pigment dispersion and aggregation characterization
- *And more....*

# Specifications & main characteristics

## SPECIFICATIONS

Particle size range (Hydrodynamic size)	0,5 nm to 10 µm (sample dependant)
Sample concentration range	Min. 0.0001% - Max. 40% (w/w) (Sample and Vasco model dependant)
Sample volume	Min.15 µL (one drop) – Max. 400µ L
Temperature control range inside the cell	10°C to 70°C +/-0,1°C
Sample Type	Aqueous & organic solvents – pH: 1-14
Accuracy / Repeatability	Better than +/- 5% with reference latex material at 100 nm
Measurement time	30 seconds to few minutes (sample dependant)
Sample Cell	Innovative and patented Dual Thickness Controller concept : extended concentration range; no need for disposable cuvette – All solvent proof ( aqueous and organic)

## SIGNAL PROCESSING

Measurement technology	Enhanced Dynamics light Scattering (DLS) – Patented Dual Thickness Controller (DTC)
Laser source	Highly reliable 65 mW (programmable) diode laser (DPSS) @635 nm coupled to optical fiber with maintain of polarization. Laser life time > 10,000 hours
Scattering angle detection	Single angle for DLS at 135° (See VASCO FLEX series for other angle)
Data processing algorithm	Proprietary inversion algorithm “PADE LAPLACE” and “Cumulants” included. Wizard Intelligent Simulation (WIS) mode for comparison with experimental results. Results display in real time in Intensity-Volume-Numbers distribution
Detectors type	Avalanche Photodiode (APD)
Correlator	Linear Type - Min. sample time 125 ns, max. delay time 100µs. Max. 1000 channels

## OPTIONS


Alternative wavelength laser	Highly reliable 30 mW diode @532 nm coupled to optical fiber with maintain of polarization
Narrow band fluorescence filter	@635 nm or @532 nm - Improves signal ratio for samples that fluoresce
High temperature control	Extends upper temperature range to 90°C
21 CFR part 11 software option	Enables an operating mode that assist with ER/ES compliance
Zeta potential measurement	See WALLIS analyzer

## GENERAL

Computer interface	USB 2.0 – Windows XP, Seven
Weight	12 kg
Dimensions	30 x 33 x 28 cm <sup>3</sup> (LxWxH)
Power supply	100-115/220-240 VAC, 50/60 Hz, 100 W max

## SYSTEM COMPLIANCE

CE certification	CE marked product - Class I laser product, EN 60825-1:2001, CDRH
ISO norm	ISO 13321 : 1996 – Particle size analysis - photon correlation spectroscopy ISO 22412 : 2008 – Particle size analysis – Dynamic Light Scattering (DLS)
Patent granted	Reference: FR2947339 (A1) – Improved particle size analyzer



Contact:  
sales@cordouan-tech.com

11, avenue de Canteranne  
33600 Pessac – France  
Tel +33 (0)5 56 15 80 45

